

Instrument color	Initial drill		Distance in mm	
	Ø 1.5 mm	Ø 2.0 mm	Tooth — implant drill	Implant drill — implant drill
	Implant Ø			
Grey	3.4-3.5	4.0	3.0	5.0
Green	3.75	4.25	3.1	5.3
Yellow	4.0-4.1-4.15	4.6	3.3	5.6
Orange	4.3	4.8	3.4	5.8
Blue	4.8-4.85	5.35	3.7	6.3

LM Implant Misura MR

by Style Italiano

For more details scan
the QR codes or go to
[www.lm-dental.com/
lm-implant-misura](http://www.lm-dental.com/lm-implant-misura)

